


passion to lead
EPSILON

EP Corsa & Astro
HF Mobile Digital Radiography

Your
Vision
in to the
Future



KW : 2.5 / 3.5 / 6 / 15 / 32

Certified & approved By:



Technical Specifications

SPECIFICATIONS	HIGH FREQUENCY MOBILE X-RAY GENERATOR						
MODEL	EP CORSA 2.5 M	EP CORSA 3.5 M	EP CORSA 4 SM	EP CORSA 6 SM	EP CORSA 6 M ASTRO 6 M	EP CORSA 15 M	EP CORSA 30 M
Max Output (KW)	2.5	3.5	4.2	6	6	15	32
kV Range	40-110	40-110	40-120	40-110	40-125	40-125	40-125
mA Max.	60	80	100	150	150	300	400
X-Ray Tube Head	Monoblock Version						
X-Ray Tube	Stationary Anode				Rotating Anode		
Digital Display	KV, mA & mAs					KV & mAs	
Control Panel	Feather Touch (Micro Processor Controller Based)						
APR Facility	Available						
Collimator	Manually / Motorized Light Beam Diaphragm (LBD)						
Power Required	1Ø 190-240 VAC 50 Hz with 15 Amps Power Socket (Required Proper Line Earthing)						
Tube Stands	Spring Balance Mobile / Counter Balance Mobile						
Battery/Capacitor Bank	Optional						

Due to continuous R & D for product development and improvement design, specifications, appearance & Data indicated might be vary from the product illustrated in the catalogue.

SPECIFICATIONS	DR FPD SERIES
Scintillator	GOS (Gadolinium Oxysulfide) / CsI (Cesium iodide)
Detector Size	14 x 17 inch / 17 x 17 inch
Type	Wired / Wireless



Epsilon Healthcare Solutions (P) Ltd.

Corporate Office : E.D Tower, F-101, J-17 D, Ramesh Marg, C-Scheme, Jaipur-302001, Raj. (INDIA) Tel : +91-141-2223773

Registered Office : 103-104, Amar Vijay Complex, Hotel Man Singh Lane, S. C. Road, Jaipur-302001, Raj. (INDIA), Tel : +91-141-2362854

Works : G-1-57, RIICO Industrial Area Extension, Phase-2, Bagru, Jaipur, Raj. (INDIA)
E-mail : marketing@epsilonhealthcare.com, Web : www.epsilonhealthcare.com

Marketed by:

Toll Free : 1800 1234 122