

passion to lead

**EPSILON**

EP Corsa HF Portable Solution



- Easy to carry
- Real compact Portable X-ray Machine
- Suitable for Home Care/ Veterinary/Medical Camp uses

**KW : 2.4**

Certified & approved By:



# Technical Specifications

## HFG for Portable / Mobile X-Ray System



APPLICATION	PORTABLE / MOBILE
MODEL	EP CORSA 2.4 P
Power Rating	2.4 kW
kV Range	40-100 kV (in 1 KV Step)
mA Range	16-35 mA 5 Step Program
mAs Range	0.4 - 100 mAs / 35 Steps
Display	kVP & mAs
X-Ray Tube	Stationary Anode
Focal Spot	1.2 x 1.2 (mm)
APR Facility	Available & Re-programmable
Line Power	1Ø 200V / 220V / 240V
Tube Stand	Optional
Battery	Optional
Weight	13 Kg (Approx)

Due to continuous R & D for product development and improvement design, specifications, appearance & Data indicated might be vary from the product illustrated in the catalogue.

## Epsilon Healthcare Solutions (P) Ltd.

**Corp. Office** : E.D Tower, F-101, J-17 D, Ramesh Marg, C-Scheme, Jaipur-302001, Raj. (INDIA) Tel : +91-141-2223773

**Regd. Office** : 103-104, Amar Vijay Complex, Hotel Man Singh Lane, S. C. Road, Jaipur-302001, Raj. (INDIA), Tel : +91-141-2362854

**Works** : G-1-57, RIICO Industrial Area Extension, Phase-2, Bagru, Jaipur-303007 Raj. (INDIA)  
E-mail : marketing@epsilonhealthcare.com, Web : www.epsilonhealthcare.com

Marketed by:

**Toll Free : 1800 1234 122**