

E3

Ventilator

- V-A/C, P-A/C, V-SIMV, P-SIMV, CPAP/PSV Ventilation mode
- Compact design
- 10.4" touch TFT LCD screen
- 3 waveforms display
- Dual CPU control
- 72 hours trend
- Internal battery support up to 90 min running time
- Optional Duo level and NIV ventilation mode
- Optional EtCO₂ monitoring



mindray
healthcare within reach

E3

Ventilator

Technical Specifications

Patient type Adults and children (bodyweight of at least 4 kg/8.8 lbs.)

Ventilation settings

Ventilation mode Assist/Control (VCV or PCV)
SIMV, SIMV/Psupp (VCV or PCV)
PSV
CPAP
DuoLevel
NIV – Mask Ventilation

Ventilation frequency (f) 1 to 100 /min (A/C), 1 to 60 /min (SIMV)

Inspiration time (Tinsp) 0.2 to 10 s

Tidal volume (VT) 0.1 to 2.0 L (Adult)/0.04 to 0.3 L (Pediatric)

Inspiratory flow 6 to 100 L/min (Adult)
6 to 30 L/min (Pediatric)

Inspiratory pressure 5 to 100 cmH₂O

PEEP Off, 1 to 45 cmH₂O

△Psupp 0 to 100 cmH₂O

Rise time for inspiratory pressure 0 to 2 s

O₂ concentration 21 to 100 Vol. %

Trigger (Pressure/Flow) 0.5 to 15 L/min / -10 to -0.5 cmH₂O

Monitored Parameters

Airway pressure Peak pressure, plateau pressure, mean airway pressure,
PEEP, min. pressure (-20 to +120 cmH₂O)

Minute volume (MV) MV, MVspont, MVleak (0 to 100 L/min)

Tidal volume (VT), (BTPS*) Inspired VT, expired VT, VTPS (0 to 4000 ml)

Breathing frequency (f) ftot, fspn, fmand. (0 to 200/bpm)

O₂ concentration (FiO₂) Inspired O₂ concentration (15 to 100 Vol. %)

Lung mechanics Resistance (cm H₂O/L/s)
Compliance (0.0 to 300 ml/cm H₂O)
RSBI (0 to 9999 l/(minx • L)
WOB (0 to 100)(J/min)
NIF (-45 to 0)(cm H₂O)
P 0.1 (-20 to 0)(cmH₂O)
PEEPi (0 to 120)(cm H₂O)

Waveforms Airway pressure-time, flow-time, volume-time,
Loops Paw-V, V-Flow, Flow-Paw
Trends 72 hours
Capnography (etCO₂) (optional)
0 to 99 mmHg with waveform

Alarms

Airway pressure High / Low
Expired minute volume High / Low
Tidal volume High
Apnea alarm time 5 to 60 s
Spontaneous breath frequency High
Inspired O₂ concentration High / Low
EtCO₂ (optional) High / Low

Performance data

Max. flow for pressure support and spontaneous breathing 120 L/min (adult)
Control principle Time cycled, volume constant, pressure-controlled
Leakage and hose system compensation automatic
Outlet for pneumatic nebulizer

Operating data

Mains power connection 100 to 240 V, 50/60 Hz
Power consumption Typically approx. 195 W
Gas supply operating pressure O₂, air: 2.7 to 6..5 bar / 39 to 94 PS/280 to 650 kPa
Physical Specifications
Dimensions ventilator (H x W x D) 330 x 390 x 365 mm (include display)
Diagonal screen size 10.4" TFT LCD color touch screen
Weight basic unit Approx. 16 kg

DISTRIBUTOR:



MINDRAY is a trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Specifications subject to changes without prior notice.
© 2011 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved.
P/N: ENG-E3-210285-20110811

mindray

Mindray is listed on the NYSE under the symbol "MR"
Mindray Building, Keji 12th Road South, High-tech Industrial Park,
Nanshan, Shenzhen 518057, P.R. China
Tel: +86 755 26582888 Fax: +86 755 26582680
E-mail: intl-market@mindray.com Website: www.mindray.com